

6.0. Aesthetic Line



ABSOLUTE Clear Aesthetic Sensation

Absolute the best aesthetic solution. When the doctor needs good low friction and aesthetic solution, the best economical aesthetic solution.

Mono-Crystalline Synthetic Sapphire Bracket

CRYSTAL CLEAR: ABSOLUTE brackets are crystal clear, and they will remain stain-free.

MAXIMUM STRENGTH: The "ABSOLUTE" has a crystal structure treated with advanced heat technology and has a torque strength three times as resistant as that of other poly-ceramic products.

SUPERIOR BOND STRENGTH AND EASIER TO DEBONDING: Absolute brackets provide the maximum bond strength and are easier to de-bond without damaging enamel.


COMFORT: The ABSOLUTE bracket is designed in a low-profile and therefore enhances patient comfort and reduces the foreign body sensation.

100% MONO-CRYSTALLINE BASE TECHNIQUE*

For more aesthetic clear and bond strength without base coating technique SLOT SCRIBE LINE MARK* for more accurate positioning.

FEATURES

- Single Crystal Structure
- 3 times better torque strength than poly ceramic products
- Easy debonding
- Small Tie-Slot Radius
- Angled Occlusal Base
- Twin Design
- Stain Resistant



6.1. M.B.T.

Maxillary	Torque	Angle	H	W	Reference Code		Reference Code	
					Slot 0,18	Slot 0,22	Right	Left
Centrals	+17°	+4°		3.6	Right	Left	Right	Left
Laterals	+10°	+8°		3.0			360-501	360-502
Cuspid w/Hook	-7°	+8°		3.5			360-503	360-504
1°-2° Bicuspid w/Hook	-7°	0°		3.5			360-505	360-506
							360-507	360-508

Mandibular	Torque	Angle	H	W	Reference Code		Reference Code	
					Slot 0,18	Slot 0,22	Right	Left
Anterior	-6°	0°		2.7	Right	Left	Right	Left
Cuspid w/Hook	-6°	+3°		3.5			360-510	360-510
1° Bicuspid w/Hook	-12°	+2°		3.5			360-511	360-512
2° Bicuspid w/Hook	-17°	+2°		3.5			360-515	360-516
							360-519	360-520

AE961-507 - Kit Patient .022 3-4-5 Hook Upr-Lwr Pkg. 20

AE961-508 - Kit Patient .022 Cuspid Hook 5-5 Upr-Lwr Pkg. 20

6.2. ROTH

Maxillary	Torque	Angle	H	W	Reference Code		Reference Code	
					Slot 0,18	Slot 0,22	Right	Left
Centrals	+12°	+5°		3.6	Right	Left	Right	Left
Laterals	+8°	+9°		3.0	360-561	360-562	360-601	360-602
Cuspid w/Hook	0°	+11°		3.5	360-563	360-564	360-603	360-604
Bicuspid w/Hook	-7°	0°		3.5	360-565	360-566	360-605	360-606
Bicuspid Universal	-7°	0°		3.5	360-567	360-568	360-607	360-608
					360-569	360-569	360-609	360-609

Mandibular	Torque	Angle	H	W	Reference Code		Reference Code	
					Slot 0,18	Slot 0,22	Right	Left
Anterior	0°	0°		2.7	Right	Left	Right	Left
Cuspid w/Hook	-11°	+5°		3.5	360-570	360-570	360-610	360-610
1° Bicuspid w/Hook	-17°	0°		3.5	360-571	360-572	360-611	360-612
2° Bicuspid w/Hook	-22°	0°		2.7	360-575	360-576	360-615	360-616
					360-579	360-580	360-617	360-618

AE961-512 - Kit Patient .018 3-4-5 w/h Upr-Lwr Pkg. 20

AE961-513 - Kit Patient .018 Cuspid Hook Upr-Lwr Pkg. 20

AE961-517 - Kit Patient .022 3-4-5 Hook Upr-Lwr Pkg. 20

AE961-518 - Kit Patient .022 Cuspid Hook 5-5 Upr-Lwr Pkg. 20