

## 6.3. TRANSLUX CERAMIC BRACKETS

**Translux is the perfect alternative for aesthetic treatment.**

**The Translux** brackets are manufactured with ceramic material of high quality. Very easy to remove due to mechanical retention in the base. **The Translux** slots are specially polished with diamonds wheels to improve the sliding properties.

**The Translux** size is the result of the design optimization, associating a reproduced dimensions part with aesthetic effect and comfort to the patient.

### ACTUAL DESIGN

With external contour and rounded edges bringing comfort to the patient and aesthetic visual.

### HIGH QUALITY BASE

With composed shaft providing a perfect adaptation.

### TRANSLUCENCY

Providing an excellent aesthetic.


### AMPLE SPACE ON THE FINS

### TORQUE IN THE BASE

### MECHANICAL RETENTION

### MAIN FEATURES

Manufactured in pure polycrystalline alumina (99,9%) and uniformity in the particles size, provide the maximum translucency, superior mechanical resistance, merging naturally on any enamel tone, once your translucency is absorbed instead reflected.



## 6.4. Translux Roth Prescription

| Maxillary       | Torque | Angle | H | W   | Reference Code |         | Reference Code |         |
|-----------------|--------|-------|---|-----|----------------|---------|----------------|---------|
|                 |        |       |   |     | Slot 0,18      |         | Slot 0,22      |         |
|                 |        |       |   |     | Right          | Left    | Right          | Left    |
| Centrals        | +12°   | 5°    |   | 3,9 | 360-401        | 360-402 | 360-451        | 360-452 |
| Laterals        | +8°    | 9°    |   | 3,3 | 360-403        | 360-404 | 360-453        | 360-454 |
| Cuspid w/Hook   | -2°    | 11°   |   | 3,3 | 360-405        | 360-406 | 360-455        | 360-456 |
| Bicuspid w/Hook | -7°    | 0°    |   | 3,3 | 360-407        | 360-408 | 360-457        | 360-458 |

| Mandibular         | Torque | Angle | H | W   | Reference Code |         | Reference Code |         |
|--------------------|--------|-------|---|-----|----------------|---------|----------------|---------|
|                    |        |       |   |     | Slot 0,18      |         | Slot 0,22      |         |
|                    |        |       |   |     | Right          | Left    | Right          | Left    |
| Anterior           | 0°     | 0°    |   | 2,8 | 360-410        | 360-410 | 360-460        | 360-460 |
| Cuspid w/Hook      | -11°   | 5°    |   | 3,3 | 360-411        | 360-412 | 360-461        | 360-462 |
| 1° Bicuspid w/Hook | -17°   | 0°    |   | 3,3 | 360-413        | 360-414 | 360-463        | 360-464 |
| 2° Bicuspid w/Hook | -22°   | 0°    |   | 3,3 | 360-417        | 360-418 | 360-469        | 360-470 |

**AE960-481** - Kit Patient .018 3-4-5 Hook 5-5 Upr-Lwr Pkg. 20

**AE960-483** - Kit Patient .022 3-4-5 Hook 5-5 Upr-Lwr Pkg. 20

**AE960-485** - Kit Patient .018 Cuspid Hook 5-5Upr-Lwr Pkg.20

**AE960-487** - Kit Patient .022 Cuspid Hook 5-5 Upr-Lwr Pkg. 20

## 6.5. Translux M.B.T. Prescription

| Maxillary       | Torque | Angle | H | W   | Reference Code |      | Reference Code |         |
|-----------------|--------|-------|---|-----|----------------|------|----------------|---------|
|                 |        |       |   |     | Slot 0,18      |      | Slot 0,22      |         |
|                 |        |       |   |     | Right          | Left | Right          | Left    |
| Centrals        | +17°   | +4°   |   | 3,9 |                |      | 360-101        | 360-102 |
| Laterals        | +10°   | +8°   |   | 3,3 |                |      | 360-103        | 360-104 |
| Cuspid w/Hook   | -7°    | +8°   |   | 3,3 |                |      | 360-105        | 360-106 |
| Bicuspid w/Hook | -7°    | 0°    |   | 3,3 |                |      | 360-107        | 360-108 |

| Mandibular         | Torque | Angle | H | W   | Reference Code |      | Reference Code |         |
|--------------------|--------|-------|---|-----|----------------|------|----------------|---------|
|                    |        |       |   |     | Slot 0,18      |      | Slot 0,22      |         |
|                    |        |       |   |     | Right          | Left | Right          | Left    |
| Anterior           | -6°    | 0°    |   | 2,8 |                |      | 360-110        | 360-110 |
| Cuspid w/Hook      | -6°    | +3°   |   | 3,3 |                |      | 360-111        | 360-112 |
| 1° Bicuspid w/Hook | -12°   | +0°   |   | 3,3 |                |      | 360-113        | 360-114 |
| 2° Bicuspid w/Hook | -17°   | +0°   |   | 3,3 |                |      | 360-117        | 360-118 |

**AE960-559** - Kit Patient .022 3-4-5 Hook 5-5 Upr. Only Pkg. 10  
**AE960-560** - Kit Patient .033 Cuspid Hook 5-5 Upr. Only Pkg.10

**AE960-561** - Kit Patient .022 3-4-5 Hook Upr-Lwr Pkg. 20  
**AE960-562** - Kit Patient .022 Cuspid Hook 5-5 Upr-Lwr Pkg. 20

## 6.6. Translux EDGEWISE STANDARD

| Maxillary | Torque | Angle | H | W   | Reference Code |         | Reference Code |         |
|-----------|--------|-------|---|-----|----------------|---------|----------------|---------|
|           |        |       |   |     | Slot 0,18      |         | Slot 0,22      |         |
|           |        |       |   |     | Right          | Left    | Right          | Left    |
| Centrals  | 0°     | 0°    |   | 3.5 | 360-001        | 360-001 | 360-021        | 360-021 |
| Laterals  | 0°     | 0°    |   | 3.5 | 360-002        | 360-002 | 360-022        | 360-022 |
| Cuspid    | 0°     | 0°    |   | 3.5 | 360-003        | 360-003 | 360-023        | 360-023 |
| Bicuspids | 0°     | 0°    |   | 3.5 | 360-005        | 360-005 | 360-025        | 360-025 |

| Mandibular                  | Torque | Angle | H | W   | Reference Code |         | Reference Code |         |
|-----------------------------|--------|-------|---|-----|----------------|---------|----------------|---------|
|                             |        |       |   |     | Slot 0,18      |         | Slot 0,22      |         |
|                             |        |       |   |     | Right          | Left    | Right          | Left    |
| Incisors                    | 0°     | 0°    |   | 2.5 | 360-010        | 360-010 | 360-020        | 360-020 |
| Cuspid                      | 0°     | 0°    |   | 3.5 | 360-003        | 360-003 | 360-023        | 360-023 |
| 1°/2° Bicuspid Univ. w/Hook | 0°     | 0°    |   | 3.5 | 360-005        | 360-005 | 360-025        | 360-025 |

**AE960-550** Kit Patient Edg. Std. .018 5-5 Upper Pkg. 20    **AE960-552** Kit Patient Edg. Std. .022 5-5 Upper Pkg. 20